

# Eppendorf Reference® 2



**1. Single-button operation:** Quick and easy tip ejection immediately following dispensing saves unergonomic side movement of thumb and actively reduces aerosols re-entering the pipette

**2. Robust design:** Sturdy upper handle including stainless steel and bumper system guarantees long service life and increased robustness

**3. Volume Setting:** Quick and secure volume setting including volume lock to prevent unwanted change. 4-digit display for accurate volume setting

**4. User adjustment:** Option to temporarily adjust the pipetting volume by up to  $\pm 2\%$  to offset inaccuracies when working with difficult liquids (e.g., warm, cold, volatile or high density liquids). Return to factory adjustment without calibration

**5. Easy to clean surface:** Smooth surface without grooves where bacteria could hide allows easy and safe cleaning

**6. Eppendorf TrackIT:** Embedded RFID data chip for heightened traceability in combination with Eppendorf TrackIT, see page 71 for details

**7. Serial number traceability:** Serial number is printed on upper and lower part to avoid a mixup of volume determining parts

**8. Fully autoclavable** without disassembly

**9. Spring-loaded tip cone:** Improves ergonomics by lowering tip attachment forces and reproducibility by ensuring consistent tip attachment

( )

T 043-235-2300 F 043-235-2301

: [www.joongsoo.co.kr](http://www.joongsoo.co.kr)



reddot design award  
winner 2013

